
NOTHING IS SIMPLE IN ONCOLOGY. NOTHING BUT THIS. idylla

Idylla™ A revolutionary, fully automated system that makes molecular testing convenient and fast. Suitable for any lab.



WE AIM TO PROVIDE RAPID ACCESS TO HIGH-LEVEL MOLECULAR TESTING ANYTIME, ANYWHERE.

THE IDYLLA™ SOLUTION - MOLECULAR TESTING SUITABLE FOR VIRTUALLY ANY LAB.



IDYLLA™ FITS ANY SIZE LAB



IDYLLA™ IS A COST-EFFECTIVE IN-HOUSE SOLUTION EVEN FOR SMALL SAMPLE VOLUMES

- Easy to implement and maintain
- No molecular infrastructure needed
- No batching needed
- No external controls needed





IDYLLA™ PROVIDES AN ULTRA-RAPID ACTIONABLE SOLUTION COMPLEMENTING COMPREHENSIVE NGS PROGRAMS

- Short turnaround time and minimal hands-on time
- Minimal sample input requirements
- On demand testing
- Seamless integration into virtually any laboratory workflow



THE REVOLUTIONARY IDYLLA™ WORKFLOW






IDYLLA™ COMPARED TO OTHER TECHNOLOGIES

MINIMAL HANDS-ON AND ASSAY TURNAROUND TIMES



REDUCED NUMBER OF INSTRUMENTS AND CONSUMABLES NEEDED



Based on workflow exercise in real-life laboratory setting. Reference technologies used: Illumina MiSeq and Qiagen Therascreen.

THE IDYLLA™ ASSAYS. MINIMAL SAMPLE REQUIREMENTS AND VERSATILE SAMPLE TYPES.



IDYLLA™ TISSUE SPECIMEN REQUIREMENTS

KRAS, NRAS & BRAF	 1 x 5 μm FFPE tissue section 1 x 10 μm FFPE tissue section Neoplastic cell content If < 10%, macrodissection needed 	50 - 600 mm² 25 - 300 mm² > 10%
BRAF	 1 x 5 μm FFPE tissue section 1 x 10 μm FFPE tissue section Neoplastic cell content If < 50%, macrodissection needed 	50 - 600 mm² 25 - 300 mm² > 50%
EGFR	 1 x 5 µm FFPE tissue section Neoplastic cell content If < 10%, macrodissection needed 	≥ 1 mm² ≥ 10%
GeneFusion	 1 x 5 µm FFPE tissue section 3 x 5 µm FFPE tissue section Neoplastic cell content If < 10%, macrodissection needed 	20 - 100 mm² < 20 mm² ≥ 10%
MSI	 1 x 5 μm FFPE tissue section 1 x 10 μm FFPE tissue section Neoplastic cell content If < 20%, macrodissection needed 	50 - 600 mm² 25 - 300 mm² > 20%



IDYLLA™ PLASMA SPECIMEN REQUIREMENT

• 1-2 ml plasma per cartridge and two validated tube types:



K²EDTA tubes,

process within 4 hours



Streck Cell-freeDNA BCTs,

process within 3 days



For Research Use Only. Not for use in diagnostic procedures.

Biocartis US Inc. 30 Montgomery Street, 9th Floor, Suite 970 Jersey City, NJ 07302 1-844-4-IDYLLA (1-844-443-9552) customerserviceUS@biocartis.com

