

# **TECHNICAL SHEET IDYLLA™ EGFR MUTATION TEST**

The **Idylla™ EGFR Mutation Test**, performed on the Biocartis Idylla<sup>™</sup> System, is an *in vitro* diagnostic test for the qualitative detection of **exon 18 (G719A/C/S), exon 21 (L858R, L861Q), exon 20 (T790M, S768I) mutations, exon 19 deletions and exon 20 insertions** in the *EGFR* oncogene. The Idylla<sup>™</sup> EGFR Mutation Test, from **sample-to-result**, starts with formalin-fixed, paraffin-embedded (FFPE) human tissue from non-small cell lung cancer (NSCLC) to liberate DNA for subsequent real-time PCR amplification and detection.

## FEATURES

EGFR mutation detection					
Exon 18	G719A G719C G719S	c.2156G>C c.2155G>T; c.2154_2155delinsTT c.2155G>A			
Exon 19	Del9	c.2238_2248delinsGC c.2239_2248delinsC c.2240_2248del cc.2239_2247del			
	Del12	c.2239_2251delinsC c.2240_2251del			
	Del15	c.2235_2249del c.2236_2250del c.2239_2253del c.2240_2254del c.2238_2252del c.2237_2251del c.2235_2252delinsAAT c.2237_2252delinsT c.2234_2248del c.2236_2253delinsCTA c.2237_2253delinsTA c.2235_2251delinsAG c.2236_2253delinsAG c.2236_2253delinsCAA c.2230_2249delinsGTCAA			



Exon 19	Del18	c.2240_2257del c.2237_2255delinsT c.2239_2256del c.2236_2253del c.2239_2258delinsCA c.2237_2254del c.2238_2255del c.2237_2257delinsTCT c.2236_2255delinsAT c.2236_2256delinsATC c.2237_2256delinsTT c.2237_2256delinsTC c.2237_2256delinsTC c.2235_2255delinsGGT			
	Del21	c.2238_2258del c.2236_2256del			
	Del24	c.2253_2276del			
Exon 20	Т790М	c.2369C>T			
	S768I	c.2303G>T			
	insG insASV9 insASV11 insSVD insH	c.2310_2311insGGT c.2308_2309insGCCAGCGTG c.2308_2311delinsCCAGCGTGGAT c.2311_2312insGCGTGGACA c.2319_2320insCAC			
Exon 21	L858R	c.2573T>G c.2573_2574delinsGT c.2573_2574delinsGA			
	L861Q	c.2582T>A			

# EGFR Total (acting as Sample Processing Control)

Specimen requirements	
Sample Type	1 x 5 μm FFPE tissue section
Neoplastic cells	$\geq$ 10%, if less macrodissection is required
Performance	
Analytical Sensitivity	LOD $\leq$ 5% for most prevalent <i>EGFR</i> mutations
	100% agreement for 10% EGFR G719S
Between Laboratory	100% agreement for 10% EGFR Del15
Reproducibility	100% agreement for 10% EGFR T790M
(600 results at 3 sites)	100% agreement for 10% EGFR L858R
	100% agreement for 10% EGFR L861Q

Between Lot Reproducibility (300 results on 3 lots) 100% agreement for 10% EGFR G719S
100% agreement for 10% EGFR Del15
100% agreement for 10% EGFR T790M
100% agreement for 10% EGFR L858R
100% agreement for 10% EGFR L861Q

Total turnaround time

Time

150 minutes

### **ACCURACY - CLINICAL PERFORMANCE EVALUATION**

J

95.9% overall agreement for the *EGFR* gene was obtained during the clinical performance evaluation comparing Idylla™ with Therascreen (Qiagen), a PCR-based reference method.

Therascreen Idylla™	DelEx19	DelEx19, T790M	L858R	L867Q	G719X	G719X, S768I	L858R, T790M	InsEx20	No mutation detected	Totals without invalids	Invalids	Totals
DelEx19	18		-	-	-	-	-	-	-	18	8	26
InsEx20	-	-	-	-	-	-	-	-	-	-	1	1
DelEx19,T790M	-	4		-	-	-	-	-	-	4	-	4
L858R	-	-	20		-	-	1	-	-	21	2	23
L861Q	-	-	-	2	-	-	-	-	-	2	-	2
G719X	-	-	-	-	1	-	-	-	-	1	-	1
G719X,S768I	-	-	-	-	-		-	-	-	1	-	1
L858R,T790M	-	-	-	-	-	-	1	-	-	1	-	1
Т790М	-	1	-	-	-	-	-	-	1	2	2	4
DelEx19, S768I	1	-	-	-	-	-	-	-	-	1	-	1
S768I	-	-	-	-	-	-	-	-	-	0	1	1
No mutation detected	-	1	2	-	1	-	-	1	66	71	30	101
Totals without invalids	19	6	22	2	2	1	2	1	67	122	-	/
Invalids	-	1	-	-	-	-	-	-	2	-	10	13
Totals	19	7	22	2	2	1	2	1	69	/	54	179

#### 95.9% overall concordance

After discordant analysis by NGS and/or ddPCR (excluding invalids) the **overall concordance** was **97.9%**.

### **IDYLLA™ EGFR PUBLICATIONS**

Check our website: https://www.biocartis.com/publications/articles-posters?cat=19 and subscribe to receive our newest publications.



Biocartis and Idylla<sup>™</sup> are registered trademarks in Europe, the United States and other countries. The Biocartis and the Idylla trademark and logo are used trademarks owned by Biocartis. Idylla<sup>™</sup> platform and Idylla<sup>™</sup> EGFR Mutation Test are CE-marked IVD's in Europe. Idylla<sup>™</sup> is available for sale in EU, USA and some other countries. Please check availability with the local Biocartis sales representative.

GL.FL009.EN.R3.07/2019