

IDYLLA™ SPECIMEN REQUIREMENTS



One of the biggest challenges in oncology biomarker testing is the ability to obtain samples of sufficient size and quality. The Idylla™ System only needs a minimal amount of sample.

SOLID BIOPSY SPECIMEN REQUIREMENTS

IVD Tests



MSI 510(k)

1 x 4 µm FFPE tissue section → 62.5-750 mm²
 1 x 5 µm FFPE tissue section → 50-600 mm²
 1 x 10 µm FFPE tissue section → 25-300 mm²
 Neoplastic cells ≥ 33% - if less, macrodissection is required

MSI-CDx (PMA)

1 x 4 µm FFPE tissue section → 62.5-750 mm²
 1 x 5 µm FFPE tissue section → 50-600 mm²
 1 x 10 µm FFPE tissue section → 25-300 mm²
 Neoplastic cells ≥ 33% - if less, macrodissection is required

RUO Assays

EGFR

1 x 5 µm FFPE tissue section → ≥ 1 mm²
 Neoplastic cells ≥ 10% - if less, macrodissection is required

GeneFusion

1 x 5 µm FFPE tissue section → ≥ 20 mm²
 3 x 5 µm FFPE tissue sections → < 20 mm²
 Neoplastic cells ≥ 10% - if less, macrodissection is required

BRAF

1 x 5 µm FFPE tissue section → 50-600 mm²
 1 x 10 µm FFPE tissue section → 25-300 mm²
 Neoplastic cells ≥ 50% - if less, macrodissection is required

KRAS NRAS-BRAF-EGFR S492R

1 x 5 µm FFPE tissue section → 50-600 mm²
 1 x 10 µm FFPE tissue section → 25-300 mm²
 Neoplastic cells ≥ 10% - if less, macrodissection is required

POLE-POLD1

1 x 5 µm FFPE tissue section → 50-600 mm²
 1 x 10 µm FFPE tissue section → 25-300 mm²
 Neoplastic cells ≥ 10% - if less, macrodissection is required

IDH1-2

50 µl extracted DNA (concentration ≥ 10 ng/µl)
 10 µl blood or bone marrow
 Maximum 3 FFPE tissue sections
 Neoplastic cells ≥ 10% - if less, macrodissection is required

PIK3CA-AKT1

1 x 5 µm FFPE tissue section → 50-600 mm²
 1 x 10 µm FFPE tissue section → 25-300 mm²
 Neoplastic cells ≥ 20% - if less, macrodissection is required

LIQUID BIOPSY SPECIMEN REQUIREMENTS¹

RUO Assays



ctEGFR

2 ml plasma collected using EDTA tubes

ctKRAS

ctNRAS-BRAF-EGFR S492R

1 ml plasma collected using either EDTA or Streck tubes

(1) K₂EDTA tubes: process within 4 hours. Streck Cell-Free DNA BCT® Tubes: process within 3 days.



Please refer to the product instructions for complete instructions on how to process samples on Idylla™.

Idylla™ Platform is listed as a class II device in the US under establishment registration 3009972873. Idylla™ MSI Test (A0160/6) is cleared in the US under K211181 for Lynch Syndrome. Idylla™ CDx MSI Test (A0220/6) is approved in the US under P250005. Idylla™ EGFR, ctEGFR, BRAF, KRAS, ctKRAS, NRAS-BRAF-EGFR S492R, ctNRAS-BRAF-EGFR S492R, POLE-POLD1 & PIK3CA-AKT1 Mutation Assays; Idylla™ GeneFusion Assay; and Idylla™ IDH1-2 Mutation Assay Kit are for Research Use Only, not for use in diagnostic procedures. For more information on the acceptable use and licenses of Idylla™ products, please visit www.biocartis.com/en/license-statements.

Biocartis and Idylla™ are registered trademarks in Europe, the US and many other countries. The Biocartis and Idylla™ trademarks and logos are used trademarks owned by Biocartis NV. Idylla™ is available for sale in Europe, the US and many other countries. Please check availability with a Biocartis representative.

© November 2025, Biocartis NV. All rights reserved.