

TECHNICAL SHEET IDYLLA™ CDx MSI TEST

For in vitro diagnostic use. For professional use only. For use on the Biocartis Idylla™ System only.

The **Idylla™ CDx MSI Test** is an in vitro diagnostic test intended for the qualitative detection of a novel panel of **seven (7) monomorphic biomarkers** for identification of microsatellite instability (MSI) in **colorectal cancer** tissue. The Idylla™ CDx MSI Test uses **formalin-fixed, paraffin embedded (FFPE) tissue** sections from human colorectal cancer (CRC), from which nucleic acids are liberated, then analyzed by PCR amplification and high-resolution melting detection. The Idylla™ CDx MSI Test automates the entire process from FFPE specimen processing to reporting of MSI status as either Microsatellite Stable (MSS) or Microsatellite Instability-High (MSI-H). The Idylla™ CDx MSI Test is intended for identification of adult patients with **MSI-H metastatic CRC**, for whom treatment with **OPDIVO® (nivolumab) in combination with YERVOY® (ipilimumab)** as a combination therapy is indicated with approved therapeutic product labeling.

FEATURES

Idylla™ CDx MSI biomarkers

ACVR2A

BTBD7

DIDO1

MRE11

RYR3

SEC31A

SULF2

Sample processing control

The MSI-specific software automatically verifies the validity of the measured fluorescence profiles. For each of the seven biomarkers, the presence of biomarker-specific amplification curves confirms successful PCR amplification and proper assay performance. These curves function as internal controls, eliminating the need for a separate external sample-processing control within the cartridge.

Specimen requirements

Sample type	1-5 FFPE tissue sections (4-10 µm)
Neoplastic cells	≥ 33%; if less, macrodissection is required
Tissue area	62.5-750 mm ² (4 µm)
	50-600 mm ² (5 µm)
	25-300 mm ² (10 µm)

Total turnaround time

Time	Approx. 150 minutes
------	---------------------

Accuracy

98.57% overall percent agreement against the OncoMate™ MSI Dx Analysis System

Idylla™ MSI Test*	OncoMate™ MSI Dx				TOTAL
	MSI-H	MSS	INVALID	NO CALL	
MSI-H	31	1***	0	3	35
MSS	1**	107	0	0	108
Invalid	0	0	0	0	0
Total	32	108	0	3	143

*Accuracy data were obtained with the Idylla™ MSI Test and are applicable to the Idylla™ CDx MSI Test.

**One (1) sample that tested MSS by Idylla™ MSI Test and MSI-H by the OncoMate™ MSI Dx Analysis System is a confirmed Lynch Syndrome case by NGS.

***One (1) sample that tested MSI-H by Idylla™ MSI Test and MSS by the OncoMate™ MSI Dx Analysis System is a confirmed Lynch Syndrome case by NGS.

Clinical performance

PFS of first line OPDIVO® (nivolumab) in combination with YERVOY® (ipilimumab) established by CheckMate 8HW

Source	Arm B: nivolumab + ipilimumab			Arm C: Chemo			Hazard Ratio***	
	N	Median** (months)	95% CI	N	Median** (months)	95% CI	Estimate	95% CI
All CTA+ Randomized	200	N.A.*	(34.30 - N.A.*)	101	6.21	(4.70 - 9.00)	0.32	(0.22 - 0.45)
Idylla™ CDx MSI Test MSI-H and CTA+	147	N.A.*	(38.44 - N.A.*)	71	6.21	(4.70 - 9.03)	0.20****	(0.12 - 0.31)
Idylla™ CDx MSI Test MSS and CTA+	30	1.81	(1.45 - 5.75)	15	7.36	(4.01 - 12.88)	1.51	(0.68 - 3.33)

*N.A.: Not available since the median PFS is not reached **Median PFS is based on Kaplan-Meier Estimates ***Hazard Ratio (HR) is Arm B over Arm C from a Cox Model stratified by tumor sidedness (left vs. right) as entered into the Interactive Response Technology (IRT) system. ****The HR is based on the Idylla™ CDx MSI Test MSI-H and CTA+ subjects. The difference between the HRs (0.20 in the table above vs. 0.21 in the drug labeling) for the first line PFS comparison of nivolumab plus ipilimumab versus chemotherapy is due to the slightly different population since the intended use population in the drug labeling is based on at least one of the two centrally confirmed tests on the CTA+ randomized subjects (Idylla™ CDx MSI Test or MMR IHC Panel pharmDx).

Catalog number

Idylla™ CDx MSI Test A0220/6

Regulatory status

Regulatory status CE-IVD – IVDR certified

Idylla™ Platform is CE-marked in Europe in compliance with EU IVD Regulation 2017/746 (IVDR) and registered in many other countries. Idylla™ CDx MSI Test (A0220/6) is CE-marked in Europe in compliance with the EU IVD Regulation 2017/746 (IVDR). The Biocartis and Idylla™ trademarks and logos are used trademarks owned by Biocartis. OPDIVO® and YERVOY® are registered trademarks of Bristol-Myers Squibb Company. Idylla™ is available for sale in Europe, the US and many other countries. Please check availability with the local Biocartis representative.

© March 2026, Biocartis NV. All rights reserved.