

FEASIBILITY STUDY OF A ctEGFR PROTOTYPE ASSAY ON THE FULLY AUTOMATED IDYLLA™ PLATFORM

BIOCARTIS

Martin Reijans¹, Sanne van Gestel¹, Elfi de Haes¹, Charlotte Vandesteene¹, Juan Manuel Castro Gomez¹, Cedric Gouedard², Stefania Patera², Samuel Murray³ and Geert Maertens¹
¹Biocartis, Generaal de Wittelaan 11B, 2800 Mechelen, Belgium;
²BioPath Innovations SA, Maroussi 15124, Athens, Greece
³BioMarker Solutions Ltd, NW52DX London, UK

INTRODUCTION

To predict the response to EGFR tyrosine kinase inhibitor (TKI) therapy in non-small cell lung cancer (NSCLC) patients, formalin-fixed paraffin-embedded (FFPE) tumor tissue is routinely tested for the presence of somatic mutations in the epidermal growth factor receptor (EGFR) gene. Sufficient tumor tissue is not always available and ctEGFR testing from plasma is an alternative approach for the detection of EGFR mutations.

METHODS

Idylla™ (Biocartis) is a fully integrated molecular detection platform that combines speed and ease of use with high sensitivity and high multiplexing capabilities. In terms of ctDNA testing, it overcomes the time-consuming step of ctDNA extraction from plasma. After the addition of proteinase K and 2 ml of human plasma into the cartridge, the complete process of ctDNA extraction, real-time PCR, data analysis and reporting is fully automated (Fig. 1). The ctEGFR prototype assay allows the detection of mutations including insertions and deletions in exons 18, 19, 20 and 21. A study was performed at BioPath Innovations SA (Athens, Greece) to compare the results obtained by the ctEGFR prototype assay with results obtained on the same plasma samples by NGS (Ion PGM platform with Ion AmpliSeq™ panel targeting EGFR exons 18, 19, 20 and 21). For each amplicon/exon, NGS target coverage was 2000 reads and the assay could detect allelic frequencies down to 2-5%. Discrepant samples were further tested with the cobas® EGFR Mutation Test v2.

RESULTS

Sixty-four NSCLC plasma samples were tested with both assays. One sample failed on NGS. Overall, 34 mutations were detected by NGS and confirmed by the ctEGFR prototype assay. In 33 samples, NGS detected no mutation. The ctEGFR prototype assay detected 8 additional mutations of which 7 mutations were detected in samples which were scored WT by NGS (Fig. 2). Retesting with the cobas EGFR Mutation Test v2 confirmed the presence of these mutations. Analytical sensitivity was assessed for 20 mutations using plasma spiked with synthetic targets. Analytical sensitivities ranging from 1 to 4% were obtained for the tested mutations (Fig. 3). Inclusivity was demonstrated for 49 mutations in total (Fig. 4). The average turnaround time of a run of the ctEGFR prototype assay was <2h 40 min and the hands-on time for loading the cartridge was <2 min.

CONCLUSION

This study shows that the Idylla™ platform enables the development of a prototype ctEGFR assay with high sensitivity and ease of use combined with a fast turnaround time for the testing of 49 relevant EGFR mutations in 2 ml of plasma.

DISCLOSURES

All authors are full-time employees of the indicated companies.

Figure 1: ctEGFR work flow




Figure 2: Concordance between EGFR genotyping of plasma with Idylla™ ctEGFR prototype assay and NGS

	NGS									
	L858R	Ex 19 del	T790M	Ex 20 insertion	S768I	L861Q	G719A/C/S	L858R/T790M	Ex19 del/T790M	No mutation
Idylla	7	12	1		1					2*
L858R										4*
Ex 19 del										1*
Ex 20 insertion										6
S768I										28
L861Q										
G719A/C/S										
L858R/T790M	1*									
Ex19 del/T790M										
No mutation										

* The presence of these mutations was confirmed by the cobas® EGFR Mutation Test v2

Figure 3: LOD table of the Idylla™ ctEGFR prototype assay

EGFR Mutation	LOD
	(% mutant in 10,000 WT cps/ml plasma)
Exon 18	
LOD c.2156G>C (G719A)	3%
LOD c.2155G>A (G719C)	1%
LOD c.2155G>T (G719S)	2%
Exon 19	
LOD c.2239_2248delinsC (del9)	2%
LOD c.2239_2251delinsC (del12)	1%
LOD c.2235_2249del (del15)	1.5%
LOD c.2236_2250del (del15)	1.5%
LOD c.2240_2257del (del18)	3%
LOD c.2238_2258del (del21)	1%
LOD c.2236-2256del (del21)	2%
LOD c.2253-2276del (del24)	4%
Exon 20	
LOD c.2369C>T (T790M)	1.5%
LOD c.2303G>T (S768I)	1%
LOD c.2307_2308insGCCAGCGTG (insASV(9))	2%
LOD c.2309_2310insCCAGCGTGGAT (insASV(11))	2%
LOD c.2311_2312insGCGTGG ACA (insSVD)	1%
LOD c.2310_2311insGGT (insG)	1%
LOD c.2319_2320insCAC (insH)	2%
Exon 21	
LOD c.2573T>G (L858R)	1%
LOD c.2582T>A (L861Q)	2%

Figure 4: EGFR mutations detected by Idylla™ ctEGFR prototype assay

EGFR exon	Mutation	AA mutation	NT mutation
Exon 18	G719A	p.Gly719Ala	c.2156G>C
	G719C	p.Gly719Cys	c.2155G>T
	G719C2	p.Gly719Cys(2)	c.2154_2155delinsTT
	G719S	p.Gly719Ser	c.2155G>A
Exon 19	Deletion 9	p.Leu747_Ala750delinsPro	c.2238_2248delinsGC
		p.Leu747_Ala750delinsSer	c.2240_2248del
		p.Leu747_Glu749del	c.2239_2247del
	Deletion 12	p.Leu747_Thr751delinsPro	c.2239_2251delinsC
		p.Leu747_Thr751delinsSer	c.2240_2251del
		p.Glu746_Ala750del	c.2235_2249del
			c.2236_2250del
			c.2239_2253del
		p.Leu747_Thr751del	c.2240_2254del
			c.2238_2252del
		p.Glu746_Thr751delinsAla	c.2237_2251del
	Deletion 15	p.Glu746_Thr751delinsIle	c.2235_2252delinsAAT
	p.Glu746_Thr751delinsVal	c.2237_2252delinsT	
	p.Lys745_Ala750delinsThr	c.2234_2248del	
	p.Glu746_Thr751delinsLeu	c.2236_2253delinsCTA	
	p.Glu746_Thr751delinsVal	c.2237_2253delinsTA	
	p.Glu746_Thr751delinsAla	c.2235_2251delinsAG	
	p.Glu746_Thr751delinsGln	c.2236_2253delinsCAA	
	p.Ile744_Ala750delinsValLys	c.2230_2249delinsGTCAA	
	p.Leu747_Pro753delinsSer	c.2240_2257del	
	p.Glu746_Ser752delinsVal	c.2237_2255delinsT	
	p.Leu747_Ser752del	c.2239_2256del	
	p.Glu746_Thr751del	c.2236_2253del	
	p.Leu747_Pro753delinsGln	c.2239_2258delinsCA	
Deletion 18	p.Glu746_Ser752delinsAla	c.2237_2254del	
	p.Glu746_Ser752delinsAsp	c.2238_2255del	
	p.Glu746_Ser752delinsIle	c.2236_2256delinsATC	
		c.2237_2256delinsTT	
	p.Glu746_Ser752delinsVal	c.2237_2256delinsTC	
		c.2235_2255delinsGGT	
Deletion 21	p.Leu747_Pro753del	c.2238_2258del	
	p.Glu746_Ser752del	c.2236_2256del	
Deletion 24	p.Ser752_Ile759del	c.2253_2276del	
T790M	p.Thr790Met	c.2369C>T	
S768I	p.Ser768Ile	c.2303G>T	
InsG	p.Asp770_Asn771insGly	c.2310_2311insGGT	
InsASV(9)	p.Val769_Asp770insAlaSerVal	c.2307_2308insGCCAGCGTG	
InsASV(11)	p.Val769_Asp770insAlaSerVal	c.2309_2310delinsCCAGCGTGGAT	
InsSVD	p.Asp770_Asn771insSerValAsp	c.2311_2312insGCGTGGACA	
InsH	p.His773_Val774insHis	c.2319_2320insCAC	
Exon 21	L858R	p.Leu858Arg	c.2573T>G
	L861Q	p.Leu861Gln	c.2573_2574delinsGA
			c.2582T>A

