

IDYLLA™ SYSTEM SPECIFICATIONS

Physical dimensions Console

w: 25.4 cm / 10 inch
h: 50.6 cm / 19.9 inch
d: 25.4 cm / 10 inch
Weight: 5.2 kg / 11.46 lb

Physical dimensions Instrument

w: 30.5 cm / 12 inch
h: 19 cm / 7.50 inch
d: 50.5 cm / 19.9 inch
Weight: 18.6 kg / 41.00 lb

Power requirements Console

Power Supply: 100-240 V; 50/60 Hz
Max. Power Consumption: 80 W
Noise Level: Max. 49 dB(A)

Power requirements Instrument

Power Supply: 100-240 V; 50/60 Hz
Max. Power Consumption: 200 W
Noise Level: Max. 54 dB(A)

Operating conditions

Temperature: 15 ... 30 °C
59 ... 86 °F
Humidity: 10 ... 80 %RH
Pressure: 80 ... 106 kPa

Cartridges

Single-use disposable cartridges
ABS plastic outer case construction



Biocartis US Inc

30 Montgomery Street, 9th Floor, Suite 970
Jersey City, NJ 07302
1-844-4-IDYLLA (1-844-443-9552)
customerserviceUS@biocartis.com