

BECAUSE
TIME
MATTERS

IDYLLA™ ctEGFR MUTATION ASSAY



Idylla™ ctEGFR Mutation Assay is for research use only. Not for use in diagnostic procedures.

FAST AND EASY ctEGFR TESTING



COMPREHENSIVE ASSAY

- Covers 49 mutations including insertions and deletions in exon 18, 19, 20 and 21 in the *EGFR* gene



EASE OF USE

- Fully automated sample-to-result process
- Walk-away system
- All PCR reagents integrated in a single cartridge
- Contamination-controlled design



FAST RESULTS

- ± 2 minutes hands-on time
- Short turnaround time ± 160 minutes from sample to result



SUITABLE FOR ANY LAB

- No special expertise required
- Minimal infrastructure needed



LIQUID BIOPSY

- Directly from 2 ml of plasma



ACCESSIBLE

- Access on demand
- No need for pre-processing or batching



THE IDYLLA™ ctEGFR MUTATION ASSAY SHOWS HIGH CONCORDANCE TO NGS¹

The Idylla™ ctEGFR Mutation Assay* compared to Next Generation Sequencing, shows **high sensitivity** and **ease of use** combined with **fast turnaround time** for the testing of 49 relevant EGFR mutations in 2 ml of NSCLC plasma samples.

- Thirty-four mutations were detected by NGS and confirmed by the Idylla™ ctEGFR prototype assay.
- An additional 8 mutations were detected by the Idylla™ ctEGFR prototype assay and confirmed by Cobas.

MUTATIONS DETECTED

Idylla™
42

NGS
34



(1) Reijans et al., ESMO 2019
*Prototype


COMPLEMENT YOUR TISSUE TESTING

Obtaining samples for biomarker testing of sufficient quantity and quality is known to be a major challenge in several cancer types including lung cancer:

- Tissue samples might be too small
- Tumor content might be too low
- Tumor might be hard to reach

Especially for these cases, testing plasma samples can be a rapid and cost-effective alternative and a valuable complement to tissue testing.

LIQUID BIOPSY WORKFLOW IDYLLA™ ctEGFR MUTATION ASSAY



Contact details

Biocartis US Inc.

30 Montgomery Street, 9th Floor, Suite 970

Jersey City, NJ 07302

1-844-4-IDYLLA (1-844-443-9552)

www.biocartis.com/US

customerserviceUS@biocartis.com

Follow us on   

For Research Use Only. Not for use in diagnostic procedures.

Copyright information

Biocartis and Idylla are registered trademarks in Europe, the United States and other countries.

The Biocartis and Idylla trademarks and logos are used trademarks owned by Biocartis.