

TECHNICAL SHEET IDYLLA™ ctNRAS-BRAF-EGFR S492R MUTATION ASSAY



The **Idylla™ ctNRAS-BRAF-EGFR S492R Mutation Assay**, performed on the Biocartis Idylla™ System, is a molecular assay for the semi-quantitative detection of 18 mutations in **codons 12, 13, 59, 61, 117, 146** of the **NRAS** gene, 5 mutations in **codon 600** of the BRAF gene and 2 mutations in **codon 492** of the EGFR gene.

The **Idylla™ ctNRAS-BRAF-EGFR S492R Mutation Assay** extracts and amplifies NRAS, BRAF and EGFR circulating tumor DNA (ctDNA) from 1 ml of plasma. The assay is a **sample-to-result** real-time PCR, intended for use in research applications only.

FEATURES

ctNRAS mutation detection		
	G12C	(c.34G>T)
	G12S	(c.34G>A)
Codon 12 (exon 2)	G12D	(c.35G>A)
	G12A	(c.35G>C)
	G12V	(c.35G>T)
Codon 13 (exon 2)	G13D	(c.38G>A)
	G13V	(c.38G>T)
	G13R	(c.37G>C)
Codon 59 (exon 3)	A59T	(c.175G>A)
	Q61K	(c.181C>A)
Codon 61 (exon 3)	Q61L	(c.182A>T)
	Q61R	(c.182A>G)
	Q61H	(c.183A>C; c.183A>T)
Codon 117 (exon 4)	K117N	(c.351G>C; c.351G>T)
Codon 146 (exon 4)	A146T	(c.436G>A)
	A146V	(c.437C>T)

NRAS Total (acting as Sample Processing Control)

ctBRAF mutation detection

Codon 600	BRAF V600E	(c.1799T>A; c.1799_1800delinsAA)
	BRAF V600D	(c.1799_1800delinsAC)
	BRAF V600K	(c.1798_1799delinsAA)
	BRAF V600R	(c.1798_1799delinsAG)

BRAF Total (acting as Sample Processing Control)

ctEGFR mutation detection

Codon 492	S492R	(c.1476C>A)
	S492R2	(c.1474A>C)

EGFR Total (acting as Sample Processing Control)

Specimen requirements

Sample Type 1 ml of plasma

Supported blood collection tubes

EDTA tubes Blood must be centrifuged within 4 hours after collection

Cell-free DNA BCT® tubes (Streck) Blood must be centrifuged within 3 days after collection

Total turnaround time

Time Approx. 110 minutes

Result Reporting

Report Qualitative genotype call and Cq values



Biocartis trademark and logo are trademarks belonging to Biocartis and are used and registered in Europe. Idylla is a registered trademark in the United States (US) and other countries. Idylla trademark and logo are used trademarks belonging to Biocartis. Idylla™ platform is a CE-marked IVD in Europe. The Idylla™ ctNRAS-BRAF-EGFR S492R Mutation Assay is available for Research Use Only (RUO), not for use in diagnostic procedures.

REF: ONCO_TS_CTNRAS_US_201812