

IDYLLA™ SYSTEM SPECIFICATIONS

Physical dimensions Console

w: 25.4 cm / 10 inch
h: 50.6 cm / 19.9 inch
d: 25.4 cm / 10 inch
Weight: 5.2 kg / 11.46 lb

Physical dimensions Instrument

w: 30.5 cm / 12 inch
h: 19 cm / 7.50 inch
d: 50.5 cm / 19.9 inch
Weight: 18.6 kg / 41.00 lb

Power requirements Console

Power Supply: 100-240 V; 50/60 Hz
Max. Power Consumption: 80 W
Noise Level: Max. 49 dB(A)

Power requirements Instrument

Power Supply: 100-240 V; 50/60 Hz
Max. Power Consumption: 200 W
Noise Level: Max. 54 dB(A)

Operating conditions

Temperature: 15 ... 30 °C
59 ... 86 °F
Humidity: 0% to 80% RH
Pressure: 80 ... 106 kPa

Cartridges

Single-use disposable cartridges
ABS plastic outer case construction

IDYLLA™ PRODUCT INFORMATION

DIAGNOSTIC PRODUCTS (IVD) - 6 CARTRIDGES/BOX	CATALOG#
Idylla™ BRAF Mutation Test	A0010/6
Idylla™ KRAS Mutation Test	A0020/6
Idylla™ NRAS-BRAF Mutation Test	A0030/6
Idylla™ EGFR Mutation Test	A0060/6
Idylla™ ctKRAS Mutation Test	A0080/6
Idylla™ ctNRAS-BRAF Mutation Test	A0090/6
Idylla™ MSI Test	A0100/6
Idylla™ SARS-CoV-2 Test	A1043/6
Idylla™ SARS-CoV-2/Flu/RSV Panel	A1050/6
SeptiCyte® RAPID	SB3010C
RESEARCH PRODUCTS (RUO) - 6 CARTRIDGES/BOX	CATALOG#
Idylla™ BRAF Mutation Assay	A0011/6
Idylla™ KRAS Mutation Assay	A0021/6
Idylla™ NRAS-BRAF-EGFR S492R Mutation Assay	A0031/6
Idylla™ EGFR Mutation Assay	A0061/6

RESEARCH PRODUCTS (RUO) - 6 CARTRIDGES/BOX	CATALOG#
Idylla™ GeneFusion Assay	A0121/6
Idylla™ ctBRAF Mutation Assay	A0071/6
Idylla™ ctKRAS Mutation Assay	A0081/6
Idylla™ ctNRAS-BRAF-EGFR S492R Mutation Assay	A0091/6
Idylla™ MSI Assay	A0101/6
Idylla™ ctEGFR Mutation Assay	A0111/6

PLATFORM (IVD) - 1 UNIT	CATALOG#
Idylla™ Instrument	P0010
Idylla™ Console	P1010

**Customer Support Department
(technical & product support)**

For product information,
technical applications and
trouble-shooting helpdesk
T +32 15 632 800
customersupport@biocartis.com

**Customer Service Department
(product ordering)**

For product ordering,
availability and
delivery details
T +32 15 632 888
customerservice@biocartis.com



Biocartis NV
 Generaal de Wittelaan 11B
 2800 Mechelen - Belgium
 T +32 15 632 888
www.biocartis.com - info@biocartis.com

The Idylla™ GeneFusion Assay, the Idylla™ SARS-CoV-2 Test, Idylla™ SARS-CoV-2/Flu/RSV Panel and SeptiCyte® RAPID contain SuperScript™ III Reverse Transcriptase and is provided subject to a license under patents or patent applications owned by or licensed to Life Technologies Corporation, which license is limited to the human diagnostic field and research field and specifically excludes applications in forensics (including human identity testing). The SuperScript™ III trademark is owned by Life Technologies Corporation.
 Patents US 7,700,339, 8,168,383, 8,481,279, 8,486,645, 8,232,060, 8,288,102, 8,377,642, 9,988,688, 9,523,130, 9,096,855, 10,526,661, 9,364,477, 9,539,254, 10,551,383 and pending US applications and all their respective foreign equivalents under license from Cell Signaling Technology, Inc.
 Idylla™ platform and Idylla™ BRAF, KRAS, NRAS-BRAF, EGFR, ctNRAS-BRAF & ctKRAS Mutation Tests and Idylla™ MSI Test are CE-marked IVD's in Europe. Idylla™ BRAF, ctBRAF, KRAS, ctKRAS, NRAS-BRAF-EGFR S492R, ctNRAS-BRAF-EGFR S492R, ctEGFR & EGFR Mutation Assays and Idylla™ MSI & GeneFusion Assays are available for Research Use Only (RUO), not for use in diagnostic procedures. The Idylla™ SARS-CoV-2 Test and Idylla™ SARS-CoV-2/Flu/RSV Panel are CE-marked IVD's in Europe, validated by Biocartis and filed for Emergency Use Authorization with the FDA. The FDA's independent review of this validation is pending. SeptiCyte® RAPID is a CE-marked IVD, developed by Immunexpress Inc and Biocartis has the exclusive distribution rights for EU and many other countries. SeptiCyte® RAPID is a registered trademark owned by Immunexpress. Idylla™ is available for sale in EU, USA and many other countries. Please check availability with the local Biocartis representative. Biocartis and Idylla are registered trademarks in Europe, the United States and other countries. The Biocartis and Idylla trademarks and logos are used trademarks owned by Biocartis. © Biocartis, September 2021.